

Heated Mitts & Liners

The therapeutic heat of the mitts soothes and nourishes the skin, nails, and cuticles while replenishing the hand's natural moisture balance. The heated mitts have a long cord between the mitts for comfort during use. Heat relieves pain caused by arthritic or other conditions as well.

Usage

1. Preheat mitten
2. Smooth hands with a softening cream, paraffin or any heated activated softening cream and massage into nails and cuticles
3. Cover hands in mitten liners (mitten with individual finger slots) or wrap in plastic wrap
4. Place hands inside preheated mitten for approximately 10 –15 minutes

Important Safety Instructions

- Read all instructions carefully.
- Do not use while sleeping.
- Do not use on an infant.
- This pad is not to be used on or by an invalid, sleeping or unconscious person or a person with poor blood circulation unless carefully attended.
- Do not use on areas of insensitive skin.
- Burns can occur regardless of control setting, check skin under pad frequently.
- Never use heated mitten without a liner of some kind.
- When using mitts, be sure to use the removable plastic liner with each unit as well as the disposable liner.