



Skin Inspection Unit

The Skin Inspection Unit allows both the technician and client to view imperfections, skin texture, and surface debris. The black light (cold quartz bulb) makes the Skin Inspection Unit a valuable tool, allowing the technician to analyze the client's skin and decide which product and procedures to use during treatment.

Parts/Accessories



Internal View



Operations/Precautions

FOR PROFESSIONAL USE ONLY

To reduce the risk of electric shock:

- Always unplug the unit when not in use.
- Do not place (or drop) into water or any other liquid.
- Do not reach for a unit that has fallen into water. Unplug immediately.

To reduce the risk of burns, electric shock, fire, or injury:

- Do not leave the Skin Inspection unit unattended when in use.
- Do not allow children to operate the Skin Inspection unit.
- Extreme caution should be used when Skin Inspection unit is used on or near invalids.
- The Skin Inspection unit should be used only for its intended purpose as described in this manual.
- Do not operate the Skin Inspection unit if:
 - It has a damaged cord or plug.
 - If it is not working properly.
 - If it has been dropped or damaged.
 - If it has been dropped into water.
- If damaged, return the Skin Inspection unit to an authorized service center for examination and repair.
- Keep the cord of the Skin Inspection unit away from heated surfaces.

- ❑ Never block the air vents of the Skin Inspection unit. Do not place it on a soft surface such as a bed or couch where air vents may be blocked. Keep the air vents free of lint, hair, etc.
- ❑ Do not use the Skin Inspection unit outdoors.
- ❑ Connect the Skin Inspection unit to a properly grounded outlet (preferably a GFCI outlet).

Instructions for Use

Setting Up Your Unit

1. Position the Skin Inspection unit on a stable, flat surface. You will want the back of the unit to be in a position where clients can sit comfortably with the curtain over their face and shoulders.
2. Release the lock located on the bottom of the Skin Inspection unit and open the lid to fully release the curtain.
3. Locate the round knob located on the side of the Skin Inspection unit. Raise the lid until it is fully extended and turn that knob counterclockwise until it locks in place.

Prior to Inspection

1. Cleanse the face and eye area or suggest that the client come in without any moisturizer or makeup.
2. Have the client sit with the curtain completely enclosing his/her head and shoulders. (The curtain has Velcro closures on the back that can be closed to keep out the light.)
3. For best results, darken the room.

Examining the Client

1. Have the client look down into the mirror.
2. After the client is in the correct position, they should close their eyes.
3. Turn the Skin Inspection unit “On” using the red On/Off Switch on the front of the unit. You will hear a fan. This is a cooling fan to remove excess heat.
4. Look through the viewing window in the front of the unit. You will now be able to see the client’s skin under ultra-violet light.
5. If client’s skin condition indicates a bright yellow film or patch, he/she may have a viral or fungal infection on their skin. Discontinue the session and refer your client to a physician or dermatologist.
6. Use the following chart to interpret the findings. Be sure to fill out a client profile card (cards are available through Universal Companies).

Color

White fluorescent
White spots
Blue-white
Purple fluorescent
Light violet
Yellow spots
Pink/Reddish
Brown
Yellow fluorescent Film

Conditions Indicated

Thick corneum layer
Dead skin cells
Healthy skin
Thin and dehydrated skin
Dehydrated
Clogged pores, comedones
Oily skin
Hyper-pigmentation
Fungal or viral infection.
(Refer client to a physician_

7. When the examination is complete, have the client open their eyes and look into the mirror briefly to point out areas of concern.
8. When the inspection is complete, turn the unit “Off” by pressing the On/Off switch.
9. When the Skin Inspection unit is not in use, unplug the unit.
10. Release the brace and fold the hood into the interior of the unit when not in use.
11. Fasten the lock and store in a secure place.

Authorized Technical Service Center:

Universal Companies, Inc.
18260 Oak Park Drive
Abingdon, VA 24210
Phone: 1-800-558-5571
Fax: 1-800-237-7199
E-Mail: info@universalcompanies.com