

1005 Cool Magnifying Lamp

【1005 Cool Magnifying Lamp manual】

This machine mainly use the magnifying function of the magnifier and cold lamp, use them to check the various small transformations of the skin surface, can timely find out the different conditions, through careful treatment and maintaining so as to keep the skin slick, delicate and stretchable.

With the unique design and magnifying function, it is also can be used in the beauty items like tattoo, grafting the eyelashes and so on. The magnifier can be used together with the bracket, or it can be used together with other instruments.

Novel design, pretty appearance. Its bracket has the protector to prevent nipping hands during operation. The magnifier adopt the lens with the usual 3 diopter magnifying rate and partially 5 diopter magnifying rate together. The soft lamplight will not make the user feel dizzy.

【Usage Methods】

- (1) Insert the electrical plug.
- (2) Turn on the power switch, turn the magnifying lamp towards the skin surface which is to be checked or operated on.
- (3) Adjusting the distance between the magnifying lamp and skin, and achieve the most clear condition.
- (4) Check the skin or start the beauty operation.

【Usage Cautions】

- (1) Forbid to use the instrument in the over humid environment, like bathroom, etc. Clean the skin surface before the beauty operation.
- (2) Turn off the power after the operation is finished.
- (3) Please draw out the plug if you won't use it for a long time. Wipe the lamp cleanly and cover it then place it inside the corner.
- (4) Please use the soft lens paper to clean the dust on the lens. Don't use other cloth or paper to clean it, otherwise the lens will be damaged.
- (5) Don't disconnect the instrument by yourself, please

ask the professionals to test or repair it if any mistake occurs.



