

Enzymatic Soaking Detergent

This product is safe for all instruments and scopes (flexible and rigid) and extends their usable lifetime. It has a neutral pH, is non-abrasive, free rinsing, and 100% biodegradable. It will not harm any metals, plastic, rubber, or corrugated tubing.

Ingredients

The Enzymatic Soaking Detergent contains a mixture of enzymes buffers and non-ionic detergents Protease, Amylase, Lipase, Carbohydrase and proprietary enzymes.

Procedure

1. Add 1 capful to soaking tray, or add contents of 16 oz. bottle to 1 gallon of water.
 2. Place treatment tips and instruments in soaking tray and allow to soak for a minimum of 15 minutes.
 3. Remove treatment tips and instruments from soaking tray and rinse.
 4. Add treatment tips to additional sanitizing solution or autoclave.
-