

MegaPeel EX[®]

Gold Series

INSTRUCTIONAL MANUAL

CONTENTS

TECHNICAL DESCRIPTION OF THE MEGAPEEL EX	4
MEGAPEEL EX PARTS AND SUPPLIES	5
INTRODUCTION	9
Chapter 1	
ASSEMBLING THE MEGAPEEL EX	11
Chapter 2	
PERFORMING A MICRODERMABRASION TREATMENT	20
MICRODERMABRASION TREATMENT GUIDELINE	22
Chapter 3	
INCORPORATING CHEMICAL PEELS WITH YOUR MICRODERMABRASION PRECISION PEEL® TREATMENT	29
Chapter 4	
INCORPORATING CHEMICAL PEELS WITH YOUR MICRODERMABRASION HYDRA BRIGHT FACIAL TREATMENT	31
Chapter 5	
TROUBLESHOOTING	34
MEGAPEEL EX® MICRODERMABRASION TREATMENT – CONSENT FORM	39
SKINCARE QUESTIONNAIRE	40
MEGAPEEL EX® TREATMENT NOTES	41
DM SKINCARE ROUTINE	43
BEFORE AND AFTER PHOTOS	48

Thank you for deciding to purchase the MegaPeel EX® Gold Series microdermabrasion system from DermaMed Solutions, LLC. We realize that it is of the utmost importance that the quality of our products and services to you match the high standards of quality that you set for the rest of your business.

This manual will guide you in using the DermaMed system. We are never more than a phone call away for any questions.

DermaMed Solutions, LLC
610-358-4447
www.DermaMedSolutions.com

MegaPeel EX[®]

Gold Series



TECHNICAL DESCRIPTION

Made in U.S.A.

Patent # 6,432,113

115 V. » 60 Hz 5A Fuse 8 Amp Type T

230 V. » 50 Hz 2.5A Fuse IEC 6.3 Amp

Class 1 Microdermabrasion Type B Applied part

This device conforms to the requirements of EN60601.1.2 for electromagnetic compatibility and will not cause or be effected by electrical noise.

Transport and storage conditions: Max50C 80 percent relative humidity

This device must be serviced only by an authorized agent. If it is serviced by an unauthorized person, the warranty will be void.

MegaPeel EX® Parts and Supplies

CRYSTALS AND FILTERS

MegaPeel EX Backup Filter

Sold in a box with 2 sets of each

Always dispose and replace crystals & waste filters simultaneously. If you do not change your waste filter along with your crystals an overflow can occur. Be sure to depress rings before disconnecting your crystal & waste filter.



MegaPeel EX Gold Crystals and Filters (CFG1000)

Sold in a box with 2 sets of each

Always dispose and replace crystals & waste filters simultaneously. If you do not waste filter along with your crystals an overflow can occur. Be sure to depress rings before disconnecting your crystal & waste filter.



DISPOSABLE TIPS, HAND-PIECES AND TUBING

Blue Tip (BT1000)

Does not come with system but available for purchase. Sold in a pack of 12

DermaMed Solutions, LLC. advises MegaPeel EX technicians to dispose of treatment tips after each client to ensure the client's safety. Before each treatment, make sure you inspect the tips for scratches or cracks. NEVER AUTOCLAVE YOUR TIPS!



Divided Clear Tip (CT1000)

Sold in a pack of 12

DermaMed Solutions, LLC. advises MegaPeel EX technicians to dispose of treatment tips after each client to ensure the client's safety. Before each treatment, make sure you inspect the tips for scratches or cracks. NEVER AUTOCLAVE YOUR TIPS!



Gold Tip (GT1000)

Sold in a pack of 12

DermaMed Solutions, LLC. advises MegaPeel EX technicians to dispose of treatment tips after each client to ensure the client's safety. Before each treatment, make sure you inspect the tips for scratches or cracks. NEVER AUTOCLAVE YOUR TIPS!



Crystal Free Tip (RCFT1000)

Does not come with system but available for purchase. Sold in a pack of 6

DermaMed Solutions, LLC. advises MegaPeel EX technicians to dispose of treatment tips after each client to ensure the client's safety. Before each treatment, make sure you inspect the tips for scratches or cracks. NEVER AUTOCLAVE YOUR TIPS!



DISPOSABLE TIPS, HAND-PIECES AND TUBING

Silver Hand-piece (HPS1000)

Most frequently used hand-piece

The clear tip with the silver hand-piece is our most frequently used combination. When using the silver hand-piece, the minimum KPA is 45. For tougher or less sensitive skin, increase the vacuum. To be replaced after 300 treatment peels.



Gold Hand-piece (HPG1000)

Medical hand-piece for aggressive treatments

The gold tip with the gold hand-piece is ideal for body and more aggressive treatments later in the patient's series. A kPa of 35 is very aggressive. The gold tip is the only tip to be used with the gold hand-piece. To be replaced after 300 treatment peels.



Tubing and Coupling (TC1000)

Recommended to be changed every 3 – 9 months depending on usage


The pink and blue tubes connect to the hand-piece. The blue tube is connected to the crystals and the pink tube is connected to the vacuum. It is important to change your tubing because it begins to wear overtime. You will know it is time to change your tubing when it becomes cloudy or the ends begin to yellow and keeps falling out of hand piece.



Spare Parts Kit (SPK1000)

For troubleshooting your MegaPeel EX

To be used in case of problems with your MegaPeel EX. Call the technical support line at 610.358.4447 for instructions. The kit of spare parts includes: 6 red O-rings, 1 clear U-shaped tube, 1 fuse, 1 pin, 1 short blue tube.



DISPOSABLE MARKETING MATERIALS

MegaPeel EX Informational Brochure (DMS2002)

Also available in Spanish & French

Promotional marketing material used to educate patients and clients about microdermabrasion. Sold individually



Microdermabrasion Gift Certificates (DMS6002)

When your customers can't decide on the "perfect gift" don't give up the sale, be ready to offer them the option of gift certificates. These Microdermabrasion Certificates are everything you need to supply and promote your microdermabrasion program.



MegaPeel EX Poster (DMS2006)

Effective way to advertise, inform your clients that you are now performing microdermabrasion.

Color poster is printed on a glossy paper stock measures 17" x 24"



MegaPeel EX Promotional Postcard (DM3003)


Promote your microdermabrasion business by mailing post cards to your clientel.

Blank backside allows for personalization and an area for mailing sticker containing your information.




MegaPeel EX Counter Card & Display Stand (DMS2004-D)

Displays the brochures in your waiting area, promoting your microdermabrasion business.



MegaPeel EX GoldManual (DMS2003-G)



MegaPeel EX Training DVD (DMTDVD)

This all-in-one video teaches you how to operate the MegaPeel EX, shows a step-by-step how to perform a microdermabrasion treatment, how to incorporate chemical peels into a microdermabrasion treatment and shows step-by-step and how to troubleshoot your system.



MegaPeel EX Informative Loop DVD (DMLDVD)

You already have the most desirable advertising situation in your waiting room: a captive audience with clients who already trust your reputation and expertise. So, take advantage of this opportunity to let them know about how you are offering microdermabrasion. In this looping promotional DVD for your waiting room, clients will learn all about microdermabrasion treatments and their benefits



DIAMOND HAND-PIECE, TIPS AND TUBING

Diamond Hand-Piece and Tip Set (DTCFKIT)

Offer clients who ask for crystal-free treatments

Diamond hand-piece and tips are embedded with high-quality industrial diamonds. The tips are available in a variety of microns/grits for different skin conditions. The combination of different crystal-free tips allows you to provide effective exfoliation for most skin types, deeper abrasion for thicker skin areas, less aggressive exfoliation, yet effective results for sensitive, fragile skin. Each tip can be easily removed from the hand piece to clean.



Diamond Hand-piece (DTH1000)

Additional hand-piece available for purchase

Extra hand-piece available for purchase



Diamond Hand-piece Tubing (DTT1000)

Recommended to be changed every 3 – 9 months depending on usage

The pink tube connect to the hand-piece. It is important to change your tubing because it begins to wear overtime.



Diamond Tubing Connector (DTC1000)

Connector used to connect Diamond Hand-piece

Connector can be replaced if necessary. Call the technical support line at 610.358.4447 for instructions.



Diamond Tip Cleaning Kit (EX-511)

For cleaning the Diamond Tips and Hand-piece

To be used to clean the Diamond Hand-pieces and Tips. Call the technical support line at 610.358.4447 for instructions.



Cavicide Disinfectant Spray (EX-903SCD)

A convenient, ready-to-use, intermediate-level surface disinfectant for the Diamond Hand-piece and Diamond Tips.



NOTE

Throughout the MegaPeel EX Manual you will notice microdermabrasion defined as mechanical exfoliation using crystals. Wherever this is stated, you may also use it to refer to microexfoliation which is done using the Diamond Hand-piece and tips.

INTRODUCTION

WHAT IS MICRODERMABRASION?

Microdermabrasion is a simple and effective form of exfoliation. It is a mechanical exfoliation that uses aluminum oxide crystals to remove the outer layer of dead or damaged skin cells, revealing a more youthful skin. Microdermabrasion is a non-invasive treatment designed to improve the overall texture, tone and clarity of the skin while restoring the skin's elasticity.

HOW DOES MICRODERMABRASION WORK?

This type of exfoliation involves gently sandblasting the skin with a high-pressure flow of crystals designed to remove the outermost layers of the epidermis (*stratum corneum*), followed by a vacuum suction that removes the crystals from the skin. While the crystals perform the exfoliation, the vacuum suction stimulates the lower layers of the skin. This repeated process stimulates the production of a new layer of skin cells, bringing with it an increased level of collagen and elastin, which further improves your skin's appearance.

HOW MANY TREATMENTS ARE NEEDED?

We recommend that people undergo a series of 4 to 6 microdermabrasion treatments. The results are cumulative. Microdermabrasion is a progressive treatment. Although 1 microdermabrasion treatment yields fresher, smoother looking skin, generally results improve substantially after the 4th, 5th or 6th treatment. The number of necessary treatments depends on the condition of the skin and the needs of the client. Initially you can treat patients every 7 to 14 days. With further treatments, they can maintain the desired results.

WHAT TYPE OF SKIN CAN BENEFIT?

All skin types and skin colors can benefit from microdermabrasion. The process rejuvenates and refreshes aging, flaky, dry, thickened, blemished, dull and sun-damaged skin. Blotchiness, uneven skin tones and texture improve. Microdermabrasion reduces the appearance of fine lines, acne scars, hyperpigmentation, chicken pox scars, brown/age spots, fine-lined wrinkles (rhytides), hyperkeratosis, large pores and other skin imperfections.

WHEN SHOULD ONE USE MICRODERMABRASION?

Use MegaPeel EX Microdermabrasion to treat:

- ◆ Face, neck and décolletage.
- ◆ Shoulders and back.
- ◆ Hands.
- ◆ Elbows and knees.

WHAT SKIN CONDITIONS RESPOND FAVORABLY TO MICRODERMABRASION?

- ◆ Acne/acne scarring.
- ◆ Dehydrated/dull skin.
- ◆ Fine lines/aging skin.
- ◆ Photo-aged skin/hyperpigmentation.
- ◆ Post inflammatory hyperpigmentation.
- ◆ Rosacea.

WHEN SHOULD ONE NOT USE MICRODERMABRASION?

Do not use microdermabrasion on people who have:

- ◆ Active acne.
- ◆ Cancer, including skin cancer.
- ◆ Eczema.
- ◆ Herpes simplex (cold sores) unless they take preventative medication before and after treatment
- ◆ Keloid scarring.
- ◆ Prescribed topical medications, such as Differin, Retin-A, Retin-A Micro, Tazorac
- ◆ Raised moles, warts, skin tags.
- ◆ Skin sensitivity.
- ◆ Sunburn.
- ◆ Uncontrolled diabetes.
- ◆ Untreated diabetes.
- ◆ Vitiligo.

Do not use microdermabrasion in:

- ◆ Areas of telangiectasias.
- ◆ Children under age 12.
- ◆ People taking Accutane within the last 6 months.
- ◆ People under the care of a physician.
- ◆ People who overuse exfoliants.
- ◆ Pregnant or nursing women.

CHAPTER 1

ASSEMBLING THE MEGAPEEL EX®

NOTE


If you plan to operate a steamer in the same room as the MegaPeel EX, call us immediately to discuss proper care. Steamers can cause crystals to clog in the MegaPeel EX.

- ◆ Take the MegaPeel EX machine out of its shipping carton.
- ◆ Store the box and packaging materials in case you need to return them.
- ◆ Put the machine on a level surface.
- ◆ Plug one end of the electrical cord into the machine and the other into a grounded, 3-pronged electrical outlet.
- ◆ Remove the lid of the MegaPeel EX and put aside.
- ◆ Inside the machine, remove the blue tube in the pink connector by pushing in on the pink ring and pull out blue tube at the same time.



CRYSTAL AND FILTER SET-UP

- ◆ Open the white box of MegaPeel EX crystals and filters.
- ◆ Take out the smaller, heavier bottle containing crystals.
- ◆ Remove the cap on the crystal container.
- ◆ Attach the container to the metal grommet at the end of the blue tube.
- ◆ Gently put this container upside down in the smaller space in the machine, on the left side.



- ◆ Take out the larger, empty waste canister with pink and blue tubes, which collects the waste from the procedures.
- ◆ Install it with the tube-side up in the space toward the right side of the machine.
- ◆ Use pressure to insert the blue tube into the blue O-ring. Push in the tube.
- ◆ Insert the pink tube into the pink O-ring.
- ◆ Replace the lid of the MegaPeel EX machine.



REMOVING AND DISCARDING USED CRYSTALS AND FILTERS

Disconnect tubes by depressing the color-coordinated connector rings. At the same time, pull out the tubes.

NOTE

The labels on the top of the Enviro-Sealed™ HEPA Cartridge (waste reservoir/filter) instruct you to change this container when you change the clean-crystal container. If you do not do this, your warranty will be void.

ASSEMBLING TUBING

- ◆ Open the 10"-by-12" cardboard box, which holds the MegaPeel EX supplies.
- ◆ Take out the pink and blue tubing. Notice that one end of the tubing has longer metal fittings than the other.
- ◆ Insert the longer grommets into the same-color O-rings on the front of the machine.
- ◆ Push to lock them into place.



NOTE

The Platinum, Gold, Silver and Bronze MegaPeel EX machines come with a silver hand-piece. The Platinum and Gold MegaPeel EX machines come with an additional gold hand-piece.

ASSEMBLING THE HAND-PIECES

Assemble the silver and gold hand-pieces the same way.

- ◆ Connect the other side of the tubing to the hand-piece.
- ◆ Put the pink end into the larger hole.
- ◆ Put the blue end into the smaller hole.
- ◆ Twist the fitting on the tubing to seal the connections. Tug the hand-piece to make sure that it is secure.



CLEANING THE HAND-PIECES

- ◆ Ream out the hand-piece using the needle in the spare-parts kit.
- ◆ Thoroughly rinse out the crystals inside the openings of the hand-piece.
- ◆ Thoroughly wash the exterior of the hand-piece.
- ◆ After treating each patient, sterilize the hand-piece by using an autoclave or alcohol.
- ◆ Replace the hand-pieces after 300 uses.

NOTE

Change the tubing and the orange pressure-locking ring every 9 months.

USING THE DISPOSABLE TIPS

From the same carton, remove the tip that you will use for treatment.

The tips determine the *depth* of the treatment. Start with the shallowest and move deeper as treatments progress.


Use the **clear tip** for most treatments, especially the first treatment. Adjust pressure to at least 45 kPa. Use the clear tip.

- ◆ Notice 2 holes on the top of the hand-piece.
- ◆ Screw the clear tip into the hand-piece, placing the divider goes between the 2 holes.
- ◆ Leave some space showing between the divider and the hand-piece.
- ◆ **CAUTION:** Do not over-twist the tip. The gap needs to remain in place.



Use the **gold tip** for patients interested in more aggressive treatments. Start with pressure of at least 35 kPa with the gold tip, which excels for areas with tougher skin, such as the back, feet and scars.

- ◆ Attach the tubing to the gold hand-piece the same way you attached it to the silver hand-piece.
- ◆ **WARNING:** Because this tip is very powerful, start the pressure at no more than 35 kPa.



NOTE

Almost all patients treated with the gold tip will experience crusting from dried serum. Crusting, which is normal, indicates that you performed the procedure properly. The crusts clear on their own in 3 to 5 days after treatment. Tell patients that they can avoid crusting if they keep the skin moist after treatment. We recommend the use of dmSkincare's argan boost serum™

You can purchase the **blue tips** separately. They are designed for those staff who are used to working at a 45 degree angle. The blue tips work for moderate to aggressive peels. Start the pressure at no less than 45 kPa.

- ◆ Attach the blue tip, screwing it on just to the point that the thin black line on the hand-set lines up with the raised line on the blue tip.



SETTING THE VACUUM PRESSURE TO START MEGAPEEL EX TREATMENT

The vacuum pressure setting determines the *power* of the treatment. The vacuum control knob is the knob on the left-hand side. The pressure plus your chosen tip determine the *depth* of the treatment as well. The MegaPeel EX uses kPa (kilopascals) as to measure vacuum pressure.

- ◆ For the 1st treatment for each patient, test the patient's reaction and comfort level.
- ◆ Use the clear tip.
- ◆ Set the vacuum pressure up to 45 kPa.
- ◆ Increase the pressure slowly if you want a more aggressive peel.
- ◆ If you are using the gold tip for a more aggressive peel, start the pressure at 35 kPa.

SETTING THE CRYSTAL FLOW

- ◆ On the outside of the machine, lower the crystal flow dial on the right by turning it to the left, counterclockwise, as far as you can without stopping it.
- ◆ Use a minimal amount of crystals – just enough to abrade the skin – to achieve the best results.

NOTE

The machine produces a little tornado of crystals in the top of the tip. Do not fill the whole side of the tip with crystals. If you have your crystals turned up too high, you will waste crystals. You will not, however, damage the machine.

- ◆ If you have questions, please visit our website or call us at 610-358-4447.

NOTE

When you place your last crystal and filter set in the system, please call and order a new box of supplies.

Diamond Hand-piece and Diamond Tips

ASSEMBLING THE DIAMOND HAND-PIECE

When using the Diamond hand-piece and Diamond tips, make sure that no other tubing is plugged into the front of the MegaPeel EX system.

- ◆ To use the Diamond hand-piece and tips, take the single pink tube and insert into the vacuum (pink) side on the front of the MegaPeel EX. The crystal (blue) side will be left without a tube.




- ◆ Insert the other end of the tubing with the metal connector into the Diamond hand-piece, making sure it is secure.



- ◆ Attach the desired Diamond tip to the end of the Diamond hand-piece.

NOTE:

Before placing the Diamond tip to the hand-piece, make sure the black o-ring is in place on the handpiece.



Daily Maintenance of Diamond Tips

CLEANING AND STERILIZATION OF TIPS IMMEDIATELY AFTER USE

Proper cleaning of the Diamond tips is imperative to your patients' health. Cleaning must be performed immediately after every treatment. It is not necessary to clean the Diamond hand-piece after each treatment.

- ◆ Immediately after each use, remove the Diamond tip from the Diamond hand-piece and rinse under hot running water. Scrub the abrasive surface of the Diamond tip with a brush (located in the cleaning kit which comes with the Diamond Tips) to remove the protein material.
- ◆ The Diamond tip should then be placed in the CaviCide1, using either a tray or an ultrasonic unit, making sure it is completely covered and let it soak in CaviCide1 for 3 minutes at room temperature (68°F/20°C).

STANDARD STEAM STERILIZATION OF TIPS

- ◆ Standard steam sterilization is based on your autoclave's maintenance parameters.
- ◆ Every sterilizer is different, please consult with your sterilizer's operating manual for recommended sterilization times and temperatures.
- ◆ Steam sterilization should always be used if the surface of the skin may be penetrated.

NOTE

General Guidelines: 30 minutes: 250° Fahrenheit at 15psi or 121° Celsius at 1kg/CM2

COLD STEAM STERILIZATION

- ◆ Most cold sterilization solutions render instruments only after a 10 hour immersion. If the instruments need to "disinfect" only, cold sterilization is ok as disinfection will take place in only 10 minutes.
- ◆ Follow all cold sterilization/disinfection protocol as mandated by your facility. Use only a facility approved disinfectant.

NOTE

General Guidelines: Soak tips in a facility approved disinfectant (such as DermaMed Solutions' Cavicide Disinfectant Spray (Item # EX-903SCD - to order please call 610.358.4447)

Maintenance of Diamond Hand-piece

CLEANING AND STERILIZATION OF HAND-PIECE

NOTE

It is not necessary to clean the Diamond hand-piece after each use

- ◆ To clean the Diamond hand-piece, first remove the Diamond hand-piece from the tubing attachment using the small screwdriver found in the Tip Box.



- ◆ Carefully insert the end of the screwdriver into the space between the Diamond hand-piece and the metal connector attached to the pink tubing. Use the screwdriver to separate the 2 pieces.



- ◆ The Diamond hand -piece should then be placed in the Cavicide1, using either a tray or an ultrasonic unit, making sure they are completely covered and let them soak in CaviCide for 3 minutes at room temperature (68°F/20°C).

STANDARD STEAM STERILIZATION OF HAND-PIECE

- ◆ Standard steam sterilization is based on your autoclave's maintenance parameters.
- ◆ Every sterilizer is different, please consult with your sterilizer's operating manual for recommended sterilization times and temperatures.
- ◆ Steam sterilization should always be used if the surface of the skin may be penetrated.

NOTE

General Guidelines: 30 minutes: 250° Fahrenheit at 15psi or 121° Celsius at 1kg/CM2

COLD STEAM STERILIZATION

- ◆ Most cold sterilization solutions render instruments only after a 10 hour immersion. If the instruments need to “disinfect” only, cold sterilization is ok as disinfection will take place in only 10 minutes.
- ◆ Follow all cold sterilization/disinfection protocol as mandated by your facility. Use only a facility approved disinfectant.

NOTE

General Guidelines: Soak tips in a facility approved disinfectant (such as DermaMed Solutions' Cavicide Disinfectant Spray (Item # EX-903SCD - to order please call 610.358.4447)

Maintenance of Diamond Tubing

CLEANING AND STERILIZATION OF CONNECTOR ON TUBING

NOTE

It is not necessary to clean the connector on tubing. The connector can NOT be autoclaved and only placed in CaviCide1

- ◆ To clean the Diamond tubing connector, depress the connector ring located at the bottom of the tubing and at the same time, pull out the pink tube.
- ◆ The tubing connector should then be placed in the Cavicide1, using either a tray or an ultrasonic unit, making sure they are completely covered and let them soak in CaviCide for 3 minutes at room temperature (68°F/20°C).
- ◆ The tubing can also be placed in the CaviCide solution while still attached to the pink tubing.

Diamond Hand-piece Treatment Guidelines

VACUUM

- ◆ The vacuum setting determines the amount of pressure between the skin and the abrasive surface.
- ◆ Low vacuum setting, 35 kPa, applies a low force.
- ◆ High vacuum setting, 45 kPa, applies a high force between the skin and the abrasive.

NOTE

Never use vacuum around the eye area. Vacuum can cause redness and swelling in this sensitive area.

- ◆ The vacuum setting can only be read when the tip end is covered.
- ◆ Cover the end of the wand and the gauge will show the current vacuum setting.
- ◆ Turning the vacuum knob clockwise will increase the vacuum pressure.

NUMBER OF PASSES

- ◆ The wands should pass over the skin surface at least twice and the passes should be overlapping to achieve a uniform result.
- ◆ The first pass is normally vertical and the second pass is horizontal.
- ◆ The speed which you move the wand is not critical.

NOTE





Refer to Chapter 2: PERFORMING A MICRODERMABRASION TREATMENT MICRODERMABRASION TREATMENT GUIDELINE in the MegaPeel EX Instructional Manual.

Diamond Tip Guide

Bonded on to the ends of the tips are natural diamonds. The diamond tips abrasiveness is measured in microns. The finer the abrasive, the less material is removed with each pass. Each tip has a number on the base which indicates the micron grit. They vary anywhere from 75-200 microns.

- ◆ The tip with 75 micron diamond abrasive is more aggressive.
- ◆ The medium tips have a diamond abrasive which read 100 and 150 microns.
- ◆ The finer tip has a diamond abrasive of 200 microns.

The diamond abrasive will not wear or become dull over time. It is very important is to keep the abrasive surface clean and to make sure you remove all skin particles.

Diamond Tip	Micron Reading	Description for Indication
	150, 75	For the body & larger surface areas, especially for acne, knees, elbow areas, and stretch marks
	150, 75	Standard tip used for face treatments.
	200, 150, 100	For scars, fine lines and smaller areas
	200,150	For scars, the nose, spot treatments of blackheads and other imperfections.

CHAPTER 2

PERFORMING A MICRODERMABRASION TREATMENT

Microdermabrasion removes dead cells from the surface of the skin. It thickens the outer layer, or epidermis, of the skin and increases the production of new skin cells, collagen and elastin. Microdermabrasion leaves the skin smooth and clean with mild redness (erythema). Complete at least 5 peels to achieve a lasting effect. The patient who follows a proper skin-care program will enjoy long-lasting improvement in the skin.

RECOMMENDED PROCEDURE:

PEEL # 1

This peel should be a conservative peel if the skin technician does not know the patient's skin and reaction to other irritants, such as glycolics, Retin-A, etc., although these products can be incorporated into a beneficial MegaPeel EX® program. Before the treatment, ask the patient if redness is acceptable and inform the patient that redness may last for a day (the pre-consultation should address this issue). Adjust setting lower if patient is apprehensive. The setting generally found to be most successful for the 1st peel is approximately 45 – 50 kPa using the clear divided tip.

PEEL # 2

When the patient comes into the office for the 2nd treatment, find out how quickly he/she recovered from the 1st treatment. Usually, the patient will report quick healing time, and with the patient's consent, the technician can be more aggressive with this peel. Increase the kPa and adjust according to the patient's tolerance. Make as many slow, even passes over the treatment area as necessary.

PEEL # 3 +

After the 2nd peel, assess the patient's recovery time. If you were aggressive, it may be necessary to space this peel 10 – 14 days after the 2nd peel. As you become more aggressive and incorporate more aggressive home care products, allow the skin to heal for up to two weeks before proceeding with another aggressive peel. The gold tip and hand-piece may be used at this point once the patient is interested in a more aggressive treatment.

CONTINUING PEELS

Follow the above guide for additional peels. All patients should have at least 5 peels, depending on the type of problems they want to address. In some cases, 10 peels are necessary for maximum results.

MAINTENANCE TREATMENTS

To maintain results and overall skin health, you should recommend that your patients return for 1-2 maintenance treatments every 3-6 months, depending on the severity of the skin challenges that were addressed in the initial treatment sessions.

POST TREATMENT PRODUCTS

As you get more aggressive with the vacuum pressure (kPa) and number of passes, you will see more erythema. As you treat the deeper layers of the epidermis, there is additional concern for hypersensitivity with some moisturizers and other skin care products. Hydrocortisone 1% (aloe cort 1% by dmSkincare) is the recommended soothing agent combined with a very non-sensitizing sunscreen (daily eclipse™ - oil free sunscreen with SPF 30 by dmSkincare).


If the patient is very irritated, use either cold compresses (excellent for calming the skin) or post-laser calming products. Continue this for two days. Do not apply any irritants, such as Retin-A or glycolics for 2 days after any peel, and use only a 1% hydrocortisone and a good moisturizing base (green tea vita-E moisturizer™ by dmSkincare).

DermaMed Solutions, LLC will be happy to recommend the proper post-skin care.

MICRODERMABRASION TREATMENT GUIDELINE


STEP 1

Begin the treatment in the center of the forehead. Pulling the skin very taut gently glide the hand piece slowly back and forth vertically from the hairline to the eyebrow moving outward toward the temples. Each pass should overlap the last without losing contact with the skin. Occasionally re-position the hands to pull the skin taut as you move along the treatment area.



STEP 2

When half of the forehead is complete move back to the center of the forehead and repeat **Step 1** on other side.




STEP 3

Have the client turn their head to the side. When treating the temple area do not lock onto the skin (this is a thin skin area that tends to be sensitive). Perform downward passes over the skin staying clear of the eye area.



STEP 4



Starting near the ear gently glide your hand piece back and forth vertically over the skin moving inward toward the nose. Remember to stay clear of the eye area; remain outside of the orbital bone. Make sure the skin is pulled very taut in the same direction the hand piece is moving.

Please note:


You may divide this area into smaller sections if necessary.

STEP 5



Have the client turn their head to the side. When treating the temple area do not lock onto the skin (this is a thin skin area that tends to be sensitive). Perform downward passes over the skin staying clear of the eye area.


STEP 6



Starting near the ear gently glide your hand piece back and forth vertically over the skin moving inward toward the nose. Remember to stay clear of the eye area; remain outside of the orbital bone. Make sure the skin is pulled very taut in the same direction the hand piece is moving.

Please note:

You may divide this area into smaller sections if necessary.




STEP 7

When treating the chin area politely ask the client to tighten up their chin by pressing their lips together firmly. Glide the hand piece vertically over the chin locking on if possible. Keep in mind that face shapes vary from person to person. If it is difficult to move vertically refer to **Steps 16 & 17** for moving the hand piece horizontally over the chin.

Please note:

You may treat along the vermilion border.

DO NOT TREAT THE LIPS.




STEP 8

Before treating the upper lip politely ask the client to tighten up the lip by pressing their lips together firmly. Gently glide the hand piece from the center of the chin outward toward the ear.

Please note:

You may treat along the vermilion border.

DO NOT TREAT THE LIPS.



STEP 9

When half of the upper lip is complete; move back to the center of the upper lip and repeat **Step 8** on other side.


Please note:

You may treat along the vermilion border.

DO NOT TREAT THE LIPS.


STEP 10

As in **Step 1** begin at the forehead. Starting at the center of the forehead gently gliding the hand piece outward toward the temples and back again toward the center of the forehead, each pass slightly overlapping the last.



STEP 11


When half of the forehead is complete move back to the center of the forehead and repeat **Step 10** on other side.



STEP 12

Have the client turn their head to the side. When treating the temple area do not lock onto the skin (this is a thin skin area that tends to be sensitive). Perform downward passes over the skin staying clear of the eye area.





STEP 13

Starting slightly below the ear gently glide the hand piece along the jaw line toward the chin and back again moving up the side of the face. Each pass should slightly overlap the last. Stay clear of the eye area; remain outside of the orbital bone. Make sure the skin is pulled very taut in the same direction you are moving the hand piece.


Please note:

You may divide this area into smaller sections if necessary.



STEP 14

Have the client turn their head to the side. When treating the temple area do not lock onto the skin (this is a thin skin area that tends to be sensitive). Perform downward passes over the skin staying clear of the eye area.



STEP 15

Starting slightly below the ear gently glide the hand piece along the jaw line toward the chin and back again moving up the side of the face. Each pass should slightly overlap the last. Stay clear of the eye area; remain outside of the orbital bone. Make sure the skin is pulled very taut in the same direction you are moving the hand piece.

Please note:

You may divide this area into smaller sections if necessary.


STEP 16

As in **Step 7**, have the client tighten up their chin by pressing their lips together firmly. Gently glide the hand piece from the center of the chin outward toward the ear.

Please note:

You may treat along the vermillion border.

DO NOT TREAT THE LIPS.




STEP 17

When half of the chin is complete move back to the center of the chin and repeat **Step 16** on other side.

Please note:

You may treat along the vermillion border.

DO NOT TREAT THE LIPS.




STEP 18


Before treating the upper lip politely ask the client to tighten up the lip by pressing their lips together firmly. Gently glide the hand piece from the center of the chin outward toward the ear.

Please note:

You may treat along the vermillion border.

DO NOT TREAT THE LIPS.






STEP 19

When half of the upper lip is complete; move back to the center of the upper lip and repeat **Step 18** on other side.

Please note:


You may treat along the vermilion border.

DO NOT TREAT THE LIPS.



STEP 20

The nose has thicker skin than the face. Slightly turn up the kPa (vacuum pressure). Gently glide the hand piece vertically and horizontally over the entire area making several passes in each direction. Include the area between the eyebrows.



STEP 21

When treating the eyes adjust your KPA to 30 - 35 (Never go above 35 KPA when treating the eye area). This area is very fragile. It is extremely important that the skin is pulled taut and away from the eye. First, move to the side of the eye. Spread the skin taut and move vertically back and forth making a few passes Next move below the eye, pull the skin down and spread the skin taut. Ideally gliding the hand piece inward making a few passes.

Please note:

DO NOT TREAT THE EYELIDS

CHAPTER 3



INCORPORATING CHEMICAL PEELS WITH YOUR MICRODERMABRASION PRECISION PEEL® TREATMENT

Chemical peels are highly effective and non-irritating treatments. Peels incorporate several chemical components and help rejuvenate the skin. You can perform peels independently or before microdermabrasion.

Precision Peel, a dmSkincare product, exfoliates dead skin cells, promotes cell renewal, stimulates circulation and promotes the growth of healthy, plump new cells. It minimizes the appearance of fine lines, scars, acne and environmental damage. Precision Peel lightens dark spots and evens skin tones. The combination of Precision Peel and MegaPeel EX microdermabrasion will provide measurable results for most skin types.

ITEMS YOU WILL NEED:

4x4 Gauze Pad

Water

botana-gel cleanser™

Peel Prep

glycolic acid 30%

treatment neutralizer

retinyl enzyme treatment

vitA-clear serum™

daily eclipse™ - oil free sunscreen SPF 30

PROCEDURE

1. Cleansing

- ◆ Apply botana-gel cleanser™ with green tea to the entire face.
- ◆ Massage gently to cleanse and remove makeup.
- ◆ Rinse with water and pat dry.

2. Peel Prep

- ◆ Apply the Peel Prep to a 4"-by- 4" gauze pad.
- ◆ Apply to the skin in a buffing motion, avoiding the eye area.
- ◆ Leave on for 2 to 3 minutes.
- ◆ Remove with a clean, dry gauze pad, buffing vigorously until the skin no longer feels tacky

3. glycolic acid 30%

- ◆ Put glycolic acid 30 % on a 4"-by- 4" gauze pad. For patients with sensitive skin, substitute lactic acid 30 %.
- ◆ Carefully apply to the skin, avoiding the eye area.
- ◆ Leave on for 2 to 3 minutes.

4. Neutralize

- ◆ Saturate a 4"-by- 4" gauze pad with treatment neutralizer.
- ◆ Generously apply to the face, avoiding the eye area.
- ◆ Rinse with water and pat dry.

5. MegaPeel EX Microdermabrasion Procedure – Refer to Chapter 2

- ◆ Dry the skin thoroughly prior to treatment.
- ◆ Make 2 passes over the skin. Treat the entire face vertically the first time, horizontally the second time.
- ◆ Pull the skin taut – in the same direction the hand-piece is moving – as you treat in that area.
- ◆ Slowly glide the hand-piece over the skin back and forth over the area.
- ◆ Dust any particles off the skin.

6. retinyl enzyme treatment

- ◆ Apply a thin layer of retinyl enzyme to a 4"-by- 4" gauze pad. For sensitive skin, substitute hydra-repair mask.
- ◆ Wipe the product over the entire face, avoiding the eye area.
- ◆ Leave on for 5 minutes.
- ◆ Use a clean, wet 4"-by- 4" gauze pad to remove the chemical.
- ◆ Rinse.
- ◆ Repeat until the skin is no longer sticky. (Do this at least 3 times.)

7. vitA-clear serum™

- ◆ Apply 3 to 4 drops of vitA-Clear Serum™ to the face using hands, avoiding the eye area. For sensitive skin, substitute green tea vita-E moisturizer.
- ◆ Do not massage the serum into the face.
- ◆ Allow to dry for 2 minutes.

8. Finish

- ◆ Apply daily eclipse™ - oil free sunscreen barrier protection to face and décolleté.

CONTRAINDICATIONS:

- ◆ Accutane use within one or two years (assessment made by physician)
- ◆ Fitzpatrick skin types IV – VI (may create hyper-pigmentation)
- ◆ Anyone with a herpetic breakout (cold sores)
- ◆ Pregnancy or lactation
- ◆ Anyone with an infectious disease
- ◆ Open sores, suspicious lesions, or basal cell carcinoma (Always refer to a physician)
- ◆ Do not treat over semi-permanent make-up

CHAPTER 4



INCORPORATING CHEMICAL PEELS WITH YOUR MICRODERMABRASION HYDRA BRIGHT FACIAL TREATMENT

Hydra Bright Facial is a procedure that will add maximum hydration and oxygen to plump tissues and brighten appearances. It incorporates D20 Solution, a liquid infusion of botanical extracts blended with deuterium oxide (“heavy water”) for intense moisture. The oxygen cream is a skin brightening cream with multiple actives and applications, combining vitamin A and glycolic acid and accelerated by peroxide which provides water and oxygen at the cellular level. This cream acts as a delivery agent for vitamin serums and creates an immediate skin lifting and visible glow. The oxygen serum uses silk amino acids and squalane to create a satiny finish on the skin and provide a super hydration and oxygenated environment to plump tissue and fill lines and wrinkles. This combination will provide measurable results for almost any skin type.

ITEMS YOU WILL NEED:

Water

botana-gel cleanser™

Steamer (optional)

peel prep

MegaPeel EX Microdermabrasion System

D20 solution

oxygen serum

oxygen cream

daily eclipse™ - oil free sunscreen SPF 30

PROCEDURE

1. Cleansing

- ◆ Apply botana-gel cleanser™ with green tea to the entire face.
- ◆ Massage gently to cleanse and remove makeup.
- ◆ Rinse with water and pat dry.

2. Peel Prep

Option 1:

- ◆ Use a steamer to soften the sebum.
- ◆ Apply peel prep with a fan brush while under steam.

- ◆ Turn steam off.
- ◆ Aggressively buff the skin dry with open-weave gauze.
- ◆ Perform extractions, if appropriate.

NOTE

Remove the microdermabrasion system from the room before applying steam.

Option 2:

- ◆ Apply peel prep to a 4"-by- 4" gauze pad.
- ◆ Apply to the skin in a buffing motion, avoiding the eye area.
- ◆ Leave on for 2 to 3 minutes.
- ◆ Remove with a clean, dry 4"-by- 4" gauze pad in a vigorous buffing motion until the skin no longer feels tacky.
- ◆ Rinse with water.
- ◆ Aggressively buff the skin dry with open-weave gauze.

3. MegaPeel EX Microdermabrasion Procedure – Refer to Chapter 2

- ◆ Dry the skin thoroughly prior to treatment.
- ◆ Make 2 passes over the skin. Treat the entire face vertically the first time, horizontally the second time.
- ◆ Pull the skin taut – in the same direction the hand-piece is moving – as you treat in that area.
- ◆ Slowly glide the hand-piece over the skin back and forth over the area.
- ◆ Dust any particles off the skin.

4. D2O solution

- ◆ Make compresses with D2O and soft, 3"-by-3" closed-weave gauze .
- ◆ Apply to the entire face, covering eyes.
- ◆ Do not cover nostrils.
- ◆ Allow the D2O to soak into the skin for 3 to 4 minutes.
- ◆ As an alternative, if oxygen tanks are available, diffuse the D2O through a spray with oxygen onto skin.

5. oxygen serum

- ◆ While skin is still damp, apply a thin, even layer of oxygen serum to entire face, including under eyes and on lips.

6. oxygen cream Mask

- ◆ Apply an even layer of oxygen cream.
- ◆ Leave it on for 5 minutes. The patient may notice temporary tingling and redness.
- ◆ Remove the residue by gently wiping with gauze dampened with D2O.

7. Hydrate

- ◆ Reapply oxygen serum, this time allowing it to absorb into the skin.

8. Finish

- ◆ Apply daily eclipse™ - oil free sunscreen to face and décolleté.

9. Home Care Treatment

After the treatment, advise the client to begin applying 3 products at home:

- ◆ C-difference™ - vitamin C lotion.
- ◆ high potenC serum™.
- ◆ collagen renewal c'reme™ with epidermal growth factor.

CONTRAINDICATIONS

- ◆ Accutane use within one or two years (assessment made by physician)
- ◆ Fitzpatrick skin types IV – VI (may create hyper-pigmentation)
- ◆ Anyone with a herpetic breakout (cold sores)
- ◆ Pregnancy or lactation
- ◆ Anyone with an infectious disease
- ◆ Open sores, suspicious lesions, or basal cell carcinoma (Always refer to a physician)
- ◆ Do not treat over semi-permanent make-up

CHAPTER 5

Before troubleshooting your MegaPeel EX be sure that the crystal canister is not empty or the waste filter is not full.

TROUBLESHOOTING

We expect that your experiences with the MegaPeel EX product and process will always run smoothly. In case you run into problems, here's a handful of suggestions.

WHAT IF THE MACHINE DOES NOT TURN ON?

Make sure that cord in the back is plugged in and the power button on the front of the machine is on. Check the fuse.

Unplug the machine.

On the back of the machine, right of the circular vent, is the black fuse box.

Use a regular, flat-blade screwdriver to pull out the front of the black fuse box.

Again using the screw driver, pop off the red fuse box. If the fuse looks blown, replace it with a new fuse.

From the spare parts kit, take a new fuse, which is a glass-and-steel tube about 1 1/4-inch long.

Using your fingers, remove the old fuse from the MegaPeel EX.

Install the new fuse.

Throw the old fuse into regular trash.



WHAT IF CRYSTALS ARE NOT COMING OUT OF THE HAND-PIECE? #1

Remove the lid.


Lift the top of the machine.

Lift the white crystal container (the smaller of the 2 containers) from the machine.

To remove the blue tube, push down with your thumb on the small, white plastic piece located on the crystal container.

Take the short blue tube from the spare parts kit.

Attach it to the metal grommet at the end of the long, blue crystal tube.



Press on the pink O-ring to remove the pink tube attached to the waste container from the machine.

Attach the pink tube from the waste filter to the other end of the short blue tube from the spare parts kit which is already connected to the blue tube coming out of the machine.

Return to the front of the machine.

Increase the vacuum pressure by turning the knob all the way to the right.



Turn on the machine.

Let the machine run for one minute.

When the machine stops, remove the short blue tube that you just installed.

Return the short blue tube to the spare parts kit.

Reconnect the long blue tube to the crystal container

Reconnect the pink tube to the pink connector inside the machine.

Run the machine.

The MegaPeel EX should now be receiving crystals.

WHAT IF CRYSTALS ARE NOT COMING OUT OF THE HAND-PIECE? #2

Disconnect the pink and blue tubes from the front of the machine.

Push down on their O-rings.

Pull out the tubes.

Take the clear, 5-inch, U-shaped tube from the spare parts kit and connect it to the O-rings on the front of the machine. The tube will look like a horseshoe in the front of the machine.

Turn the vacuum knob all the way to the right so the vacuum is at its highest pressure.



Turn on the machine.

If you see crystals going through the clear tube, turn the machine off.

Remove the clear U-shaped tube.

Replace all connections: pink tube, blue tube, hand-piece and tip.

Turn the machine on.

Occlude the tip 5 to 10 times.

The MegaPeel EX should now be receiving crystals.

WHAT IF CRYSTALS ARE STILL NOT COMING OUT OF THE HAND-PIECE? #3

Remove the lid..

Check the pink and blue tubing which you attach the hand-piece to and see if there is a clog.

If you see a clog, replace the tubing.

If you do not see a clog, disconnect the hand piece.

Use the needle in the spare-parts kit to unclog the holes.

Knock the hand piece forcefully against a hard surface.

If you still do not see crystals in the tube, repeat the process with the short blue-tube as mentioned in Step #1 on page 25


Unplug the machine.

Lift the top of the machine.

Lift the white crystal container (the smaller of the 2 containers) from the machine.

To remove the blue tube, push down with your thumb on the small, white plastic piece located on the crystal container.

Take the short blue tube from the spare parts kit.



Attach it to the metal grommet at the end of the long, blue crystal tube.

Press on the pink O-ring to remove the pink tube attached to the waste container from the machine.



Attach the pink tube from the waste filter to the other end of the short blue tube from the spare parts kit which is already connected to the blue tube coming out of the machine.

Return to the front of the machine.

Increase the vacuum pressure by turning the knob all the way to the right.



Turn on the machine.

Let the machine run for one minute.

When the machine stops, remove the short blue tube that you just installed.

Return the short blue tube to the spare parts kit.

Reconnect the long blue tube to the crystal container

Reconnect the pink tube to the waste-filter connection.

Run the machine.

The MegaPeel EX should now be receiving crystals.

IF YOU HAVE A STEAMER IN THE ROOM:

Steam can cause clogging problems. We do not recommend keeping your machine in the same room as a steamer, if you do not have a choice the following information can help you keep your system running smoothly:

- ◆ Please do not store your crystals in the same room as a steamer.
- ◆ Keep the protective cap and recap your crystals at night.
- ◆ Store crystals in a dry place (not under the sink or in the same treatment room as the steamer).

Remove the MegaPeel EX from the room when performing steam treatments. Steam remains in the air even after the steam treatment is complete

IF STEAM STILL SEEMS TO BE AN ISSUE:

After your last treatment of the day:

- ◆ Take the short blue tube from the spare parts kit.
- ◆ Connect it to the blue crystal tube and to the pink tube on the waste filter
- ◆ Leave the blue tube connected overnight.
- ◆ Run the system for several minutes in the morning before replacing the crystal canister. This will introduce dry air into the system.

If after doing all the above clearings and crystals are not coming out of the hand-piece, call DermaMed's Technical Support line, 610.358-4447. We will help you make the machine work properly.

MegaPeel EX[®] Microdermabrasion Treatment

CONSENT FORM

Prior to receiving this treatment, I have been candid in revealing any condition that may have a bearing on this procedure, such as, pregnancy, recent facial peels or surgery, allergies, tendencies to cold sores and fever blisters, use of Retin-A, Accutane or Hormones.

I understand there may be some degree of minor discomfort, i.e., scratchiness, itchiness.

I understand there are no guarantees to this procedure.

I understand that to achieve maximum results, I will need several ongoing treatments and will need to use a daily product over a period of time.

I understand that the possibility of irritation and redness exists and that I should notify my skin care professional when irritation persists.

I will follow the home care program specifically designed for me without changing or adding any products without consulting with my skin care professional.

I have read the enclosed consultation and understand the contents.

I agree to all of the above to have this treatment performed on me and will follow all prescribed directions regarding post peel care.

PRINT NAME: _____

ADDRESS: _____

SIGNATURE: _____ DATE: _____

Skin Care Questionnaire

Date: _____
Name: _____ Birthdate: _____
Address: _____
City: _____ State: _____ Zip: _____
Home Phone: _____ Work Phone: _____
Referred by: _____

PERSONAL DATA:

Smoker: (circle one) no yes Pregnant: (circle one) **no yes**
Cosmetic surgery: (circle one) **no yes** If yes, when: _____
Define procedure(s): _____
Medication: (circle one) **no yes** If yes, what kind(s)? _____
Any health problems? (circle one) **no yes** If yes, explain: _____
Any allergic reactions to medication? (circle one) **no yes** If yes, describe: _____
Do you have any allergies? (circle one) **no yes**
Do you suntan? (circle one) **no yes**
Do you use sunscreen? (circle one) **no yes**
Please name the brand of products you are currently using:
Cleanser: _____ Toner: _____
Moisturizer: _____ Scrub: _____
Mask: _____ Buff Puff: _____
Other: _____
Have you ever used Retin-A? (circle one) **no yes** If yes, what strength? _____
Have you ever been treated with Phenol or Trichloroacetic acid? (circle one) **no yes**
Have you ever used Hydroquinone (skin lightener)? (circle one) **no yes**
Have you ever been on Accutane? (circle one) **no yes** If yes, when? _____
Have you ever had *herpes, hives, cold sores, fever blisters, keloids*? Circle all that apply
If yes, when? _____
How would you characterize your skin: (circle one) **Sensitive Rough Dry Oily/Acne-prone**
If you had one complaint about your skin, what would it be? _____
Describe your skin in three words: _____
Additional comments/concerns: _____

MegaPeel EX Treatment Notes

Today's Date: _____ Treatment #: _____ Date of Last Treatment: _____

Examined client's skin and observed the following: _____

Cleansed client's skin with: _____

Prepared the skin for microdermabrasion with: _____

Pre-peel pad or solution

Allowed skin to: (circle one) *Air dry* *Fan dry* *Cloth dry*

Determined the hand-piece to be used: (circle one) *Gold* *Silver*

Determined the tip to be used: (circle one) *Gold* *Blue* *Clear divided*

Determined the vacuum power by occluding the tip – vacuum power used: _____

Determined crystal flow by turning the crystal regulating knob ____ full revolutions down from full crystal flow

Tested vacuum power on client's mid-forehead – results: _____

i.e. Slight redness, pin-point bleeding, excoriation, scratch

If above results were more or less than desired, re-adjustment of vacuum power and/or crystal flow were performed as follows: _____

Proceeded with microdermabrasion treatment (use facial chart to record vacuum power, crystal flow, hand piece and tip used) making two passes in opposite directions on the following areas (for face): Right forehead, Left forehead, Right temple area, Left temple area, Right cheek area, Left cheek area, nose, chin.

Turned the vacuum power down to: _____, to microdermabrade the right eye area, left eye area, right side of neck, mid-neck and left side of neck.

Incorporated a chemical peel/procedure: _____

Post microdermabrasion treatment:

- ◆ Removed any excess crystal residue
- ◆ Sprayed with mineral water or wiped with damp cotton
- ◆ Cleansed with: _____ (name of product)
- ◆ Applied 1% Hydrocortisone (if skin appeared very irritated) yes no
- ◆ Moisturized with: _____ (name of product)
- ◆ Applied mineral sunscreen of at least 30 SPF: yes no
- ◆ Handed client written, post treatment protocol

Attach copy of Post Treatment Protocol

SKIN CONDITIONS:

- ___ Sun Damage
- ___ Hard bumps under skin
- ___ Dry Patches
- ___ Brown Spots
- ___ Milia: *Red* *White*
- ___ Broken Blood Vessels
- ___ Upper Lip Lines: *Deep* *Fine* Clogged Pores
- ___ Pimples: *Often* *Sometimes*
- ___ Freckles
- ___ Excessive Oiliness
- ___ Wrinkles
- ___ Acne
- ___ Blackheads

SKIN TYPE

- Normal/Combination
- Oily
- Dry

Notes and Observations



HYDRATING ROUTINE

TREATMENT SCHEDULE

(Dry, Sensitive, and Mature Skin)

Name: _____

DATE					
TIME					

RECOMMENDED PROGRAM:

- ◆ MegaPeel EX® treatment every 7 to 14 days for a series of 6-8 treatments
- ◆ Maintenance treatments every 4 to 6 weeks
- ◆ Immediately after each treatment: **daily eclipse** sunscreen

HYDRATING ROUTINE:

In the morning, cleanse the face (recommended **gentle cleansing bar**). Allow skin to remain slightly damp, then apply a pea size of **argan boost serum**, smoothing over face and neck. The serum may be applied to eye area, keeping away from lash line. Follow with **C-difference -vitamin C lotion** and sunscreen (recommended daily eclipse.)

At night, cleanse the face and apply a pea size of argan boost serum on damp skin. Follow with **green tea vita-E moisturizer**.

DO NOT SPEND TIME IN THE SUN IF YOU ARE USING RETINYL PRODUCTS

ACNE ROUTINE

TREATMENT SCHEDULE

(Acne-Prone and Oily Skin)

Name: _____

DATE					
TIME					

RECOMMENDED PROGRAM:

- ◆ MegaPeel EX® treatment every 7 to 14 days for a series of 6-8 treatments
- ◆ Maintenance treatments every 4 to 6 weeks
- ◆ Immediately after each treatment: **daily eclipse** sunscreen

CLARIFYING ROUTINE:

In the morning, cleanse the face (recommended **repair bar**). After cleansing, apply **green tea toner** with a cotton pad. Keep out of eyes. Apply sunscreen (recommended **daily eclipse**).

At night, cleanse skin (recommended **botana gel cleanser**) and pat dry. Apply a pea-size amount of **vitA-clear serum** to the face. Keep out of eye area. Allow to dry for two (2) minutes before applying a dime-size amount of **AHA salicylic solution** to the face. Keep out of eye area.

DO NOT SPEND TIME IN THE SUN IF YOU ARE USING RETINYL PRODUCTS

LIGHTENING ROUTINE

TREATMENT SCHEDULE

(Hyperpigmented and Uneven Toned Skin)

Name: _____

DATE					
TIME					

RECOMMENDED PROGRAM:

- ◆ MegaPeel EX® treatment every 7 to 14 days for a series of 6-8 treatments
- ◆ Maintenance treatments every 4 to 6 weeks
- ◆ Immediately after each treatment: **daily eclipse** sunscreen

LIGHTENING ROUTINE:

In the morning, cleanse the face (recommended **botana-gel cleanser – with green tea**). Apply a dime-size amount of **C-difference** and follow with sun protection (recommended **daily eclipse**).

At night apply pea size amount of **clari-tone serum** and allow to dry for two (2) minutes. Follow with a dime-size application of **AHA salicylic solution** to the entire area. Keep both products out of eye area.

DO NOT SPEND TIME IN THE SUN IF YOU ARE USING RETINYL PRODUCTS

ANTI-AGING ROUTINE

TREATMENT SCHEDULE

(Mature/Aged Skin and Fine Lines)

Name: _____

DATE					
TIME					

RECOMMENDED PROGRAM:

- ◆ MegaPeel EX® treatment every 7 to 14 days for a series of 6-8 treatments
- ◆ Maintenance treatments every 4 to 6 weeks
- ◆ Immediately after each treatment: **daily eclipse** sunscreen

ANTI-AGING ROUTINE:

In the morning, cleanse the face (recommended **botana-gel cleanser – with green tea**). Apply a pea size amount of **high potenC serum** to face and another to neck area. Blend into skin before applying moisture products (recommend **C-difference™ -vitamin C lotion** or **green tea vita-E moisturizer™**) and sunscreen (recommended **daily eclipse™**). Do not apply **high potenC serum™** to thin skin of upper eyelids.

At night, cleanse the face and apply dime-sized amount of **collagen renewal C'reme** for nighttime repair. Three times a week (or as directed by your skincare technician) apply a quarter-size amount of **retinyl enzyme treatment** evenly to clean face and neck and leave on for (8-10) minutes, remove with warm water and gently rub with damp cloth. Repeat removal and rinse process three times. Pat skin dry and follow with your recommended serum or moisture products. Keep out of eye area.

DO NOT SPEND TIME IN THE SUN IF YOU ARE USING RETINYL PRODUCTS

POST TREATMENT/SENSITIVE ROUTINE

TREATMENT SCHEDULE

(Post Treatment, Irritated Skin & Rosacea)

Name: _____

DATE					
TIME					

RECOMMENDED PROGRAM:


- ◆ MegaPeel EX® treatment every 7 to 14 days for a series of 6-8 treatments
- ◆ Maintenance treatments every 4 to 6 weeks
- ◆ Immediately after each treatment: **daily eclipse** sunscreen

POST-TREATMENT ROUTINE:


In the morning cleanse the face using the **gentle cleansing bar**. Create a lather with water, then apply to skin and gently cleanse with fingertips. Rinse well. May also be used in evening. Apply a small amount of **argan boost serum** to skin, adding water and smooth over face, keeping out of tear ducts and away from lash line (may be used under eyes, on lips and throat). For temporary relief of itching, discomfort, inflammation and rashes that might be due to stimulating facial products or treatments apply liberal amount of **derma renewal gel**, as needed for soothing and moisture. Keep out of eye area.

DO NOT SPEND TIME IN THE SUN IF YOU ARE USING RETINYL PRODUCTS

Before and After Photos



Before and After Photos



Before and After Photos



Before and After Photos

