bellabaci™

FACIAL CUPS



Start saying goodbye to wrinkles, aging skin, and broken capillaries with the hand-squeezed Bellabaci Facial Massage Cups.

Facial Cupping is a safe non-invasive, non-surgical, luxurious express facial in the comfort of your own home.

Add this easy and ultra-relaxing step to your morning ritual for instantly more radiant skin.

FACIAL CUPPING BENEFITS

- Increased blood flow will boost collagen formation, leading to diminished wrinkles
- · Combine with active aromatherapy oils to penetrate deeper, increasing results
- Excess fluid will drain faster, relieving under-eye bags and puffiness
- A deeper cleanse will be achieved, removing comedones
- Broken facial capillaries will vanish as fresh blood will remove the de-oxygenated debris
- Prepares the skin for extraction, thereby lessening the possibility of trauma to the skin

FACIAL CUP BENEFITS

- For the general usage on the face, neck, and decollate
- Available in Soft and Hard Strength
- · Always start with your soft cup for approx. 4 weeks before graduating to hard cup
- Sensitive or very loose skin should always only be treated by the soft cup

INDICATIONS

- Wrinkles
- Sinus problems
- Decreased elasticity
- Poor blood circulation
- Fine Lines
- Broken Capillaries/Telangiectasia
- · Sallow Skin
- Poor lymphatic drainage (under the eye etc.)
- Congestion i.e. comedones

CONTRAINDICATIONS

Active Acne

Cancer

• Thrombosis

- · Rosacea / Weeping eczema
- Botox for the 1st five days after injections Sunburn
- Open cuts or abrasions

DIRECTIONS

The Facial cups can be used in 3 different methods Firstly: Use the Facial Cups in conjunction with your Cleanser: Lather your skin well to provide ample slip to move the cup about. Use quick circular motions to deeply cleanse and remove any debris. Rinse well and repeat if necessary.

Secondly: Use your Facial Cups in conjunction with your enzyme exfoliant: Apply and make sure there is enough slip for the cup to move freely. Add

water if necessary and remember to keep the cup moving in circular motions.

Thirdly: Use your Facial Cups in conjunction with one of the facial Genies or the Nautilus Luminescence anti-aging balm. Achieve deeper product penetration by combining your cups with your choice of product. Use slow movements, such as big circles and draining from the midline to the temples, ears, or jawline, as this is where your facial lymph nodes are situated.

CARING FOR YOUR CUPS

- Wash with a disinfectant liquid soap after each use.
- Do not use any mineral oil or mineral oil containing product with your Bellabaci cups.
- Take note of seasonal changes in temperature or if you use
- air-conditioning, as too hot or too cold will affect the pliability of your cups. They will become harder in cold temperatures and softer in warmer
- To combat this, either place them in hot water to make them softer, or put them in the fridge for a few minutes to make them harder.