

**Product Name**

Super Sunnies® Protective Eyewear

**Procedure:**

1. Untie elastic and push ends through openings provided on each Super Sunnies Eyeshield eye cup, working from back to front.
  2. Tie a double knot at each end of elastic.
  3. Adjust for head size and proper fit and tie another double knot to shorten the elastic, if needed. Cut off excess elastic.
  4. Never use a Super Sunnies Eyeshield that is broken, cracked, or missing a lens, due to eye safety hazard of sunlamp products. Replace worn or broken elastic.
  5. Super Sunnies are safe for use with sunlamps and medical UV light treatment devices, provided that such devices are used in strict accordance with original manufacturer's directions.
  6. Super Sunnies are in full compliance with U.S. Federal Law 21 CFR 1040.20 performance standards for sunlamp products.
  7. Proper fit of the Super Sunnies is necessary for the safety of the sunlamp user. Adjust product by bending nose bridge or shortening elastic so that no light leaks around the edge of the Super Sunnies where it meets your face.
  8. Never use a sunlamp without eye protection.
-