Spares

Only those spare parts supplied or specified by Prestige Medical should be used in the maintenance of the autoclave. Use of unauthorised parts will invalidate any warranty given and may adversely affect the performance and safety of the unit.

Accessories

A range of accessories are available for your autoclave as described below and pictured on page 76. Contact your supplier for full details.

1	-	219294	-	Lifting	device
0		010000		0	1 -

2 - 219293 - General Instrument Tray

3 - 219292 - Standard Basket4 - 219295 - 'V' Support

5 - 229310 - Dental Instrument Tray

6 - 219291 - Cassette box

7 - 219290 - Cassette Rack8 - 219296 - Extended Basket

9 - 219500 - Green Silicone Sealing Gasket

10 - 219277 - TST Indicator strips

11 - 219258 - Cord Set UK

12 - 219297 - Cord Set EURO

13 - 219299 - Cord Set UL

14 - 219509 - Extra large basket (210030/34 only)

15 - 219282 - Fuses (2x8+2x10 amp - 230v only)

Warranty

Prestige Medical shall, in the first 12 months from the date of purchase, repair or replace free of charge any parts which prove to be defective in workmanship and/or materials. The heating element (only) is covered by a lifetime guarantee.

Prestige Medical shall not be so liable in the event that the purchaser has failed to adhere to the instructions contained herein or if the autoclave has been abused, interfered with, altered, repaired or serviced by any unauthorised party this may also result in the protection provided by the equipment being impaired.

This warranty excludes the gasket, all internal furniture and consumables.

Consumer's statutory rights are not affected.

Product decontamination.

Should the unit require repair, it must be decontaminated in accordance with a recognised procedure prior to return or on-site repair. A statement of equipment contamination status must be available with the product. (Details of a suitable procedure are available on request).

Approvals.

Note: Please ask for details of approvals specific to your model.

IEC 601-1	NEMKO/SEMKO	UL 544
ESCHLE	CE Mark - Medical Device Directive	CSA
BS3970 Part 4 (pressure vessel only)	(92/42/EEC)	UDT
DIN VDE 0750	DIN 58946 Pt. 8	SEV

Packaging.

All Prestige Medical packaging materials used are recyclable, please dispose of accordingly.

PLEASE NOTE: English is the original language for the purposes of these instructions. All other languages are translations from the English text.