PROCEDURE CHAIR 300 SERIES

Coming Early Fall 2013





300 SERIES CHAIRS Specifications

STANDARD	
Motion I	Powered 19" - 35" (48 - 89cm) Height Range
Motion 2	Powered I5°/5° Tilting
Motion 3	Powered 0-80° Backrest
Motion 4	Powered 90° Leg Rest
Hand Control	All motions and memory positions
Paper Roll Holder	Accommodates up to 18" (46cm) rolls (rolls not included)
Table Voltage Options	US or International electronics
Patient Capacity	550 lbs. (250kg)
3 Section Top	74" (188cm) long (80" (203cm) long with Head Rest extended) x 27" (69cm) wide
Backrest	26" long x 27" wide (66 x 69cm)
Seat Section	23" long x 27" wide (58 x 69cm)
Leg Rest Section	13" long x 23" wide (33 x 58cm)
Padding	3" (8cm) or 4" (10cm) Comfort Foam™
Fabric	Premium, PVC free TerraTouch™
Enclosure	Black, Grey, Beige, Cool White
Chair Weight	300 lbs. (136kg)
Safety Listings*	Will meet FDA, ETL, CE requirements

 $[\]ensuremath{^{*}}$ See pg 14 for more information

OPTIONS	
Adjustable Head Rest	11" long \times 10" wide (28 \times 25cm), adjustable angle & height, removable
Large Head Rest	11.5" long x 20" (29cm x 51cm) wide, removable
Retractable Foot Rests (Stirrups)	Pull out from under seat section, fully adjustable
Foot Control	Supports all powered motions
Procedure Tray	12.25" x 10.25" (31 x 26cm), Quarts/1.9 Liters
Swivel Base	240° of rotation (adds 1.5" (4cm) to height range)
T-rail	Available for factory mount on one or both sides of the seat
ACCESSORIES	
T-rail Adapter	Attaches to the backrest for any US format T-rail attachable accessory.
Procedure Chair Face Rest Crescent	For use when prone positioning is required.
Arm Rests	Auto level design adjusts to backrest angle, swing out 90°, 24" L x 5"W (61 x 13cm), removable

for more options & accessories see pgs 14-18