

Kozi™ Shoulder Wrap

This ultra soft shoulder wrap is designed to cover the contours of the neck and shoulders. Perfect for soothing tight, tired muscles, it will stay in place while your client walks, stands, or sits. Aromatherapy benefits are increased when it is heated in the microwave. Packaged in a sturdy attractive cardboard retail box that features a Velcro® opening and a convenient nylon handle for carrying.

Ingredients

Lavender, valerian root, white willow, chamomile, rosemary, peppermint, spearmint, hops, yellow dock, lemongrass, cinnamon, yarrow, premium flax seeds, and other natural materials.

Procedure

- ◆ Never allow the Shoulder Wrap to come in direct contact with the skin.
 - ◆ Protect the Shoulder Wrap and the body by applying a cloth to the area to be treated before applying the B Shoulder Wrap.
 - ◆ Always test the Shoulder Wrap for comfort before applying.
 - ◆ Heat Therapy – Warm in microwave in 30 second increments until desired temperature is desired. ***Do not overheat: Can scorch or catch fire!***
 - ◆ Cold Therapy. Wrap the Shoulder Wrap in a towel and place in the freezer for 1-2 hours.
 - ◆ When not using the Shoulder Wrap for prolonged periods of time, place it in a sealed container or ziplock bag to increase its longevity.
-