## **Rosa Graf Concentrated Treatment Ampoules**

Each Multipulle bottle contains 30 ml of a concentrated complex of active ingredients. Each ampoule is sufficient for 15-20 treatment. They can be mixed with each other for special treatment purposes. Each bottle comes with a 2 ml syringe to measure out the correct amount needed for a treatment. If unopened, the Multipulle concentrates can be kept for 3 years. Once opened, they should be used within 3 months.

## Features/Benefits

reatures/ Denents		
Allergy		Ideal for intensive care of hypersensitive skin.
		Has a positive effect on allergic skin reactions, acne, superficial skiin
		injuries, sunburn, and photodermatosis.
		Reliably and naturally alleviates irritations.
Cacao		Formulated to infuse rich moisture into the skin for a smoother, younger
		looking complexion.
Calming		Effective on irritated skin - e.g. excess exposure to the sun or over reaction
		to an exfoliation or heat treatment.
		Contains natural sweet almond protein which cools and panthenol which
		eases irritation.
Caviar		Stabilizes and improves skin's repair mechanism.
		Treats existing atrophies.
Clearing		Treat blemished skin which tends to form comedones, papulae, and
C		pimples.
		Should be used right after a thorough cleansing.
		Contains lactic acid which has a keratolytic effect, cleansing cucumber
		extract and urea, which promotes improved moisture balance.
Collastin		Designed for skin that is prone to atrophies.
		Stimulates skin's immune system to prevent environmental aging and
		soften existing lines.
		Helps increase the skin's capacity to store moisture.
Couperose		Reduces the appearance of couperose - fine red-blue veins just below the
1		skin's surface.
Deep Cleansing		Changes the sebum structure and softens comedones with fatty acid esters.
Effecting		Combats the appearance of fine lines and wrinkles on the face.
		Can also be used for bust and décolleté treatments.
		Helps firm and tone.
Moisture		For dry skin conditions and keratosis.
		Redresses skin's moisture balance.
Oxygenating		For tired, pale skin which tends to look dull.
78 8		Contains an effective plant-based stimulant for cell respiration which
		increases the presence of oxygen in the cells and helps revitalize them.
Q10		Emulsion based, Q10 is a unique radical trap.
		Extracted from plant cultures that are identical to that of the skin.
		Helps prevent aging due to environmental factors.
		Helps balance out minor skin damage and stabilize the skin's immune
		system.
Repair		Promotes regeneration
1		Increases skin's ability to store moisture for a smother complexion

## Usage

Multipulle Concentrates can be introduced into the skin in a variety of ways:

- Massaged into the skin apply contents of syringe directly on the skin and massage in for about 5 minutes.
- Massaged into the skin under steam empty the contents of the syringe onto your hands and apply it on the skin. Use the steamer for 5 minutes.
- Introduced using ultrasound mix the contacts of the syringe with a gel suitable for use with ultrasound in a small container, and then apply the mixture to the skin using a cosmetic brush. Proceed with the ultrasound treatment.
- Introduced using iontophoresis with a positive polarity empty the contents of the syringe into your hands, and then apply it on the skin. Proceed with the galvanic treatment with the unit set to positive polarity.

If desired, one-half of the 2 ml portion can be applied under a masque and the remainder can be massaged into the skin before applying your regular massage cream.