

DEEP FLOW ATTACHMENT

Experience next-level recovery with the Tend Flow, a smoother way to work through tension. Powered by our patented **MagnoMotion™ Technology**, the Flow Attachment delivers ultra-smooth, consistent vibration across larger muscle groups while still keeping the focused accuracy that the Tend Deep is known for.

Engineered for Bigger, Better Coverage

Made from **aerospace-grade aluminium**, the Flow features a free-moving roller-ball head that glides effortlessly across the skin. The loose ball creates a continuous rolling motion, helping the vibration travel evenly and comfortably into tight or overworked areas.

- Ideal for quads, hamstrings, glutes, calves, lats and back
- Wide coverage with precision feel
- Designed for long, flowing recovery sessions

The roller ball can be **removed and cleaned** easily with warm, soapy water, quick, hygienic, and built for daily use.

Designed to Work With Tend Deep

The Flow is an attachment only.

It connects directly to the **Tend Deep battery base**, which provides the power and vibration frequency required for optimal performance.

For anyone new to TEND, the Flow can be purchased as part of the **Tend Deep Flow Bundle**, which includes the Deep battery base and an additional attachment.

The Science Behind Tend

TEND devices use **Focal Vibration Therapy (FVT)** — a research-backed method that sends high-frequency vibrations directly into targeted areas. This helps:

- Improve circulation
- Reduce inflammation
- Calm protective muscle tension
- Support mobility and movement quality
- Speed up post-activity recovery

By stimulating local tissues and the nervous system, FVT helps the body switch out of “tight and guarded” mode and into a more relaxed, responsive state — without force, pressure, or strain.

Recover faster. Move better. Feel better with TEND.

